

Samson PowerDrive

PD-6[®]



Advanced Technologies That Move You

The adjustable PD-6 handlebars are ergonomically designed for comfort and easy operation. The PD-6 exclusive dynamic brake system takes away the fear and physical stress of steep slopes.



Samson PowerDrive PD-6[®]

Infuse your wheelchair with energy!

The electric Power Drive-6 attaches to **any manual wheelchair** in seconds and takes on the chore of pushing the chair. Powerful enough for steep inclines and long range hauls, the PD-6 is extremely lightweight, easy to handle, easy to attach and detach. It's ideal for transportation in the trunk of a car.

- Easily clamps on and detaches from wheelchair.
- Forward and backward operation.
- Long range and plenty of power with the 12 A/h battery.
- Adjustable height, ergonomically designed handlebars.
- Dynamic, electromagnetic break system.
- Battery is easily removed for recharging.

Technical Specifications:

Carrying Capacity	115 kg	253 lbs
-------------------	--------	---------

Weight:

Total Weight (7A/h 12A/h)	17 kg/19 kg	37 lbs/42 lbs
---------------------------	-------------	---------------

Weight of PD-6 without batteries	11 kg	24 lbs
----------------------------------	-------	--------

Weight of the battery 7A/h	6 kg	13 lbs
----------------------------	------	--------

Weight of the battery 12A/h	8 kg	18 lbs
-----------------------------	------	--------

Charger	100-240 VAC	100-240 VAC
---------	-------------	-------------

Tiller type	Adjustable	Adjustable
-------------	------------	------------

Speed	Adjustable up to	6 km/h	4 mile/h
--------------	------------------	--------	----------

Range:

7 Ah battery	up to 10 km	up to 6 miles
--------------	-------------	---------------

12 Ah battery	up to 15 km	up to 9 miles
---------------	-------------	---------------



PD-6



PD-6 attached to a wheelchair



PD-6 in a car trunk

Weights & measurements are approximate

